



Cordless Handheld Massager with Heat Model: WS906

USER MANUAL

EN

INTRODUCTION


To fully benefit from the massage experience offered by this product, please read the user manual before first use and keep the manual for future reference.

SAFETY PRECAUTIONS

- Individuals who are pregnant, fitted with pacemaker, suffering from diabetes, phlebitis and/or thrombosis, or is at an increased risk of developing blood clots due to recent surgery or prolonged bed rest should consult their physician before using the product.
- Do not immerse the massager in water, use the product when taking shower or bathing, or in an environment which has dripping or splashing liquid around. If you spilled liquid on the product, turn off the product and dry it immediately with a soft, lint-free cloth.
- Do not subject the unit to excessive force, shock, dust, temperature or humidity.
- Do not cover the massager with anything, e.g. a blanket or clothing, while in use; excessively heating the massager may eventually cause skin burn, fire hazard or electricity leakage.
- Do not use the product on erupted, swollen, inflamed skin areas or when you have a skin disease.
- For personal hygiene reasons, wipe the product clean with a soft and slightly damp cloth after each use.
- If you are feeling drowsy, stop using the product and do not leave the product operating when you are sleeping.
- It is not recommended to receive massage before bedtime; the massager may induce stimulating effect which may cause insomnia.
- Do not use the product after drinking alcohol or when you are on medication.
- Do not prolong the massage on a single point of your body. The product is designed to massage in motions on an area of your body.
- Keep out of reach of children. Not suitable for children's use.
- The appliance is not intended for use by persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- The massager is intended for external use only.
- Do not use the product on or near eye area, any highly sensitive areas or a body part which has a poor blood circulation.
- Check the massager regularly for possible damage. Stop using when you find the item is damaged in anyway or not functioning properly.
- Do not attempt to repair the appliance yourself. Do not tamper with the unit's internal components. This invalidates the warranty.
- Do not lift the massager by pulling the attaching power cord or use the power cord as a handling unit.
- Placement of the product on certain types of timber surface may bring damage to its finish; Oregon Scientific will not be responsible for which damage.
- Consult the maintenance section on this manual for care instructions.
- Remove dust regularly with a dry cloth. Accumulated dust poses a potential fire hazard.
- The product is intended for indoor use only.
- When disposing, ensure the product is collected for proper waste treatments.
- Images shown in this manual may differ from the actual display.
- The contents of this manual may not be reproduced without the permission of the manufacturer.

NOTE: The technical specifications for this product and the contents of the user manual are subject to change without notice.

OVERVIEW



- Heat massage surface
- Massager head
- Power / Vibration speed button:** Turn on and select the intensity of vibration
- Heat button:** Activate heat function
- Massager stand
- High vibration speed indicator light (White light)
- Power on / Low vibration speed indicator light (White light)
- Heat indicator light (Red light)
- Battery indicator light (Steady white light for charging/ Flashing white light for low-battery)
- DC jack (At the bottom of the massager)
- DC / USB plug

BATTERY CHARGING

The unit should be charged before the first use and when the battery indicator light (green) is flashing (which indicates the battery is low).

NOTE: Do not operate the massager when charging.

- You can use the stand as the charging station (as illustrated in the Overview diagram or charge the unit with the DC / USB cable alone.
- Insert the DC plug of the DC / USB cable to the slot at the bottom of the stand as illustrated below.



- Insert the USB plug of the DC / USB cable to a PC or laptop USB outlet. Alternatively, you can use the AC/DC recharging adaptor to charge the unit from a electrical socket outlet.



- Flip to open the silicon cover of the DC slot. Place the massager on the stand as illustrated in Overview. Press the unit to the stand to make sure the cable is connected properly. The battery indicator light will stay on.
TIP: Make sure the stand is sat on a flat and supportive surface when pressing the

IMPORTANT: To maintain the battery life of the massager, it is recommended only charge the unit when the battery is low. Overcharging will adversely affect performance.

USING THE MASSAGER

IMPORTANT

- This massager is not a medical or precise device, and is not intended to replace professional medical care. It is only intended for leisure massaging. If you suffer from musculoskeletal disorder, consult a qualified physician / physiotherapist.
- Do not massage on a single point of your body for too long. The product is designed to massage in motions on an area of your body.

To massage at low vibration mode:

- Press **Power / Vibration speed** once. The white indicator light will appear under the massager head.
- To turn off, press **Power / Vibration speed** twice.

To massage at high vibration mode:

- Press **Power / Vibration speed** twice. Two white lights indicators will appear.
- To turn off, press **Power / Vibration speed** again.

NOTE:

- The massager switches off every 30 minutes automatically.
- A fully-charged massager can operate approximate 3,5 hours of massaging.

HEAT FUNCTION

- To activate heat function, press **Heat button** once. A red indicator light will appear. It will take a few seconds for the heating surface to heat up.
- To turn off, press **Heat button** again.

NOTE:

- The heat function switches off every 15 minutes automatically.
- A fully-charged massager can operate approximate 4,5 hours of heating.

IMPORTANT: Misapplication of the heat may cause serious injury to your skin. Do not use on sensitive or damaged skin or on areas with poor blood circulation. Do not allow children and incapacitated persons use this function.

MAINTENANCE

- Ensure the unit is unplugged and cool before cleaning.
- For personal hygiene reasons, wipe the product clean firstly with a slightly damp cloth then with a soft dry cloth before and after each use. This unit is not water proof. Soaking it in water will cause damage to the electric components. Do not clean the unit with sharp kitchen utensils or scouring pads, or corrosive materials.
- Do not clean the massager with materials or tools that could cause damage, e.g. sharp kitchen utensils or scouring pads.
- Do not switch on or operating the product when cleaning.
- Store in a dry and cool place.
- Do not expose the product to direct sunlight or store in a damp place.
- Always switch off after use.

ABOUT OREGON SCIENTIFIC

Visit our website (www.oregonscientific.com) to learn more about Oregon Scientific products. If you're in the US and would like to contact our Customer Care department directly, please visit: www2.oregonscientific.com/service/support.asp. For international inquiries, please visit: www2.oregonscientific.com/about/international.asp.

EU-DECLARATION OF CONFORMITY

Hereby, Oregon Scientific, declares that this Cordless Handheld Massager (Model: WS906) is in compliance with EMC directive 2004/108/EC. A copy of the Declaration of Conformity is available on request via our Oregon Scientific Service.

